



## Presenting the **PrinCE-C 700** flexible CE systems

## **Smart laboratory solutions from Prince Technologies.**

At Prince Technologies we've been working hard to provide cutting edge compact instruments to the scientific community. We're proud to present the fully automated PrinCE-C 700 series, powerful integrated Capillary Electro-Chromatography systems, which allow the user more hands off time.

Flexible, modular and suitable for almost any application the PrinCE-C 700 couples to any external CE detectors, including UV/VIS, CD, LIF and or MS without the need for any special tools.

Large buffer vials eradicate the need for a replenishment system reducing down-time, making the PrinCE-C 700 systems unique in the CE field.

**Flexible** 

- Ultra short effective capillary length from 7.2cm
- Outlet end injection
- Sample/buffer cooling to a broad temperature range

Modular

- Suitable for a broad range of applications including (R&D, Food, Forensic, Pharmaceutical and Life Science)
- Choice of 2 types of diode array detectors (512 & 1024)

Practical

- Large sample injection volume
- Dual pressurisation
- Increase analysis time
- Complete temperature control of capillary
- Powerful automated analysis software package
- Unrivaled reproducibility and sample injection

## **PrinCE-C 700 Series Specifications**

	ented Dynamic Compression Injection (DCI) to generate and apply a range of positive and ed ramping. Electrokinetic injection features controlled ramping in voltage and/or current mode
Duration	Up to 650 minutes, 0.01 minutes resolution
Hydrodynamic	-180 to 250 mbar, 1mbar resolution
Electrokinetic	-30 to +30kV
Modes	Current and/or voltage
Ramping	Programmable voltage and/or current
Flushing system	
Flush range	0 to +2500 mbar, up to 10 bar by external pressure
Features	Inlet or outlet and dual pressurisation
Autosampler/fraction collecto User replaceable sample and bu position carousel	r ffer segments. Automatic randomly optimised choice of positions for inlet and outlet. 30/48
Temperature control	4-40°C
Type of vials	300µl inserts, 4ml vials with resealing snap starburst caps including Eppendorf micro centrifuge tubes
Sample volume	Minimum 10μl, maximum 4ml
Buffer volume	Maximum 4ml
Capillary compartment	
Complete temperature control	5°C below ambient up to 60°C
Capillary	Standard: minimum effective capillary length 8.5cm, minimum total length 30cm Optional: minimum effective capillary length 7.2cm, minimum total length 30cm
I <b>ntegrated detector</b> Real time UV, VIS NIR diode array	v detector
PrinCE-C 760	190-610nm (512 diodes)
PrinCE-C 770	190-1010nm (1024 diodes)
Wavelength accuracy	< 0.3nm
Repetitive accuracy	< 0.1nm
Noise level	< 3* 10-5AU
ntegration time	6ms - 10s
	ion and data analysis software package. Peak purity comparing spectrum and the end of a peak, ation database, combined peak lists different wavelengths, import and export of data and others
Channels	Up to 8 channels
Raw data	Simultaneously detects signals, pressure, monitor current, outlet electrode current, oven temperature, tray temperature and more
Unique features	Current and/or voltage controlled method step durations (for on-line pre-concentration). Positive or negative pressure during electrophoresis. Outlet end injections for samples or reagents. Automatic cross-over
Features	Programmable time parameters within each step: start time within the step pressure voltage and/or current ramping to set-point + set-point, external events and extended marker programming
Real time display	Inlet, outlet, pressure, voltage, current, methods, external events, oven and tray temperature
Real time display	Inlet, outlet, pressure, voltage, current, methods, external events, oven and tray temperature 115/230V, 50/60Hz, 300VA



Weight

45kg